

**Data sheet**

This data sheet provided for information only

**2X UHotMas<sup>CF</sup>Mix-1525****RESEARCH USE ONLY****Ready-to-use amplification MasterMix**

<b>Cat.No.</b>	<b>Pack</b>	<b>Conc.</b>
<b>MMUH-100</b>	<b>100 rnx</b>	<b>2 X</b>
<b>MMUH-500</b>	<b>500 rnx</b>	<b>2 X</b>

**DESCRIPTION**

**2X UHotMas<sup>CF</sup>Mix** is a ready-to-use premix of all components for amplification of target DNA, contains **stabilizer/ enhancer**, which improves thermostabilization of enzyme during PCR amplification and storage.

**2X UHotMas<sup>CF</sup>Mix** contains **blend of antibodies-blocked SmarTaq polymerase and special additive**, which are not active at ambient temperature (during PCR set-up) and activated automatically during the first PCR cycle at the temperature >70°C, preventing miss-priming and other artifacts formation.

It is no need for prolonged heating for activation of enzyme for PCR.

**2X UHotMas<sup>CF</sup>Mix** contains optimized buffer reagents which greatly improve specificity of PCR with **complex, low-copy number DNA templates, multiplex PCR, "real-time" PCR**, allowing to use very small initial quantities of DNA template.

One' can use an appropriate volume of **2X UHotMas<sup>CF</sup>Mix** for amplification reaction, depending on total final reaction volume.

Just place it into the tube/plate adds primers and template of choice mix all components and run PCR.

**After PCR reaction running mix 5-10µl of reaction mixture with appropriate volume of "Loading Buffer" (for non-"real-time" mode PCR), apply to the gel and run electrophoresis.**

**Stability:**

**2X UHotMas<sup>CF</sup>Mix** stable for 24 months at -20°C, or for 6 months at +4°C storage without freezing.

**CONTENT:**

**1X: UHotSmarTaq** Polymerase

**0.2mM** each of dNTP's

**1,5 mM** MgCL<sub>2</sub>

Reaction Buffer components

**Stabilizer/enhancer**

**Recommended PCR assay**

<b>50µl PCR assay</b>		<b>Final Conc.</b>
25µl	<b>2X UHotMas<sup>CF</sup>Mix 1525</b>	1X
0.2-1µM	each Primer	
Variable*	DNA Template	
To 50µl	PCR Grade Water	

\*- depending on DNA template initial concentration

**APPLICATIONS:**

- Routine PCR
- Primer extension
- Real-Time PCR (all types)
- Low-copy PCR (UHotSmarTaq Polymerase)
- Multiplex PCR

**STORAGE CONDITIONS :**

Store **2XUHotMas<sup>CF</sup>Mix** at -20°C (for long-term storage).

**SHIPPING CONDITIONS:**

Should be shipped at ambient temperature

For long distance shipments preferably in **Blue Ice**

## General Protocol for amplification with 2X UHotMas<sup>CF</sup>Mix -1525

### Add and mix the following components:

Component	50 $\mu$ L reactions	25 $\mu$ L reactions	Final concentration
PCR grade Water	Up to 50 $\mu$ L	Up to 25 $\mu$ L	
<b>2XUHotMasCFMix-1525</b>	25 $\mu$ L	12.5 $\mu$ L	1X
Primers			0.3-0.5 $\mu$ M each
Template DNA	optionally	optionally	1-10ng

In some cases, we recommend to optimize Mg concentration in the range 2.0-3.0mM  
We recommend using 25 $\mu$ L reaction for the PCR with **2X UHotMas<sup>CF</sup>Mix**

### Cycling Protocol:

Cycle step	3-step amplification		Cycles
	T $^{\circ}$ C	Time	
<b>Initial Denaturation</b>	<b>95<math>^{\circ}</math>C</b>	<b>1-2 min</b>	<b>1</b>
<b>Denaturation</b>	<b>95<math>^{\circ}</math>C</b>	<b>10 S</b>	
<b>Annealing</b>	<b>55-66*</b>	<b>5-10 S</b>	<b>25-35</b>
<b>Extension</b>	<b>72<math>^{\circ}</math>C</b>	<b>15 Sec/Kb**</b>	
<b>Final extension</b>	<b>72<math>^{\circ}</math>C</b>	<b>1-2 min</b>	<b>1</b>
	<b>4<math>^{\circ}</math>C</b>	<b>hold</b>	

\*Optimal T<sub>m</sub> for the primer pair recommended as T<sub>m</sub> of the lower primer, for the standard oligos <20nt.

To optimize amplification we recommend using gradient PCR amplification, to reach final amplification conditions in the short time.

\*\*For non-complex DNA templates (plasmid DNA, phage DNA, BAC clone) extension time could be reduced up to 15 sec/Kb.

For complex DNA, templates (human DNA) strongly recommended to apply Extension time as 30 sec/Kb for the targets more than 1,5Kb