

Data sheet

This data sheet provided for information only

HS-Fuzz DNA POLYMERASE**RESEARCH USE ONLY****Recombinant HS-Fuzz DNA polymerase**

(Deoxynucleosidetriphosphate: DNA Deoxynucleosidyltransferase E.C. 2.7.7.7.)

Cat.No.	Pack	Conc.
DHSF-100	100 Units	2 U/μl
DHSF-500	500 Units	2 U/μl

Stability:

Shelf life 18 months if store at -20°C

UNIT DEFINITION

One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTP's into acid-insoluble form in 30 minutes at 74°C under assay conditions.

ASSOCIATED ACTIVITIES:

Endonuclease and exonuclease activities were not detectable after 2 and 1 hours incubation, respectively, of 1 μg lambda DNA and 0.22 μg of EcoR I digested lambda DNA, respectively, at 72°C in the presence of 15-20 units of **UltaSmarTaq** DNA polymerase.

STORAGE AND DILUTION BUFFER:

20mM Tris-HCL (pH 8.0);100mM KCL;0.1mM EDTA; 50% glycerol,1% Tween-20, Stabilizer
OR (Glycerol Free)

20mM Tris-HCL (pH 8.0);100mM KCL;0.1mM EDTA;1% Tween-20, Stabilizer

AMPLIFICATION BUFFERS:

HS-Fuzz DNA Polymerase is supplied with **2.5x Uni Buffer-1500**, containing 3.75 mM

STORAGE CONDITIONS :

Store **HS-Fuzz** DNA Polymerase at -20°C for long term storage

SHIPPING CONDITIONS:

Should be shipped at ambient temperature. For long distance shipments preferably in **Blue Ice**

SOURCE:

E.coli strain expressing the cloned **HS-Fuzz** DNA Polymerase gene. **Anti-Pfu** mouse monoclonal antibodies

DESCRIPTION

HS-Fuzz DNA Polymerase is a unique artificial enzyme created on the basis of intellectual protein design planning by genetic engineering technique.

HS-Fuzz enzymatic activity at ambient temperature is blocked by highly specific monoclonal antibodies against **HS-Fuzz** polymerase. Antibodies allows to switch enzymatic activity during first denaturation step of PCR run, but keeps **HS-Fuzz** inactive during PCR reaction setup.

This feature results in improvement of specificity of PCR and increasing of yield of PCR product, preventing miss priming and PCR artifacts formation – primer-dimers.

The enzyme possess **high fidelity** feature. The processivity of the enzyme is very high, so the combination of processivity with fidelity results in dramatically increased yield of PCR products, very high sensitivity of PCR tests, ability to amplify “difficult” templates. This dramatic increase in processivity results not only in shorter extension times, but in more robust amplification and the ability to amplify long templates really fast.

HS-Fuzz DNA Polymerase possesses the 5'->3' DNA polymerase activity, 3'->5' exonuclease activity and temperature-dependent strand-displacement activity and generates blunt ends in the amplification products.

HS-Fuzz Provides:

- **Superior fidelity - 50x improvement compared to Taq polymerase;**
- **Excellent performance across a wide range of “difficult” templates;**
- **Long range amplification of complex targets - > 10 kb from genomic DNA;**
- **High speed - reduce reaction times.**
- **dUTP poisoning resistance**
- **Resistance to blood containing DNA samples (up to 20% of blood)**

APPLICATIONS:

- Hi-Fi PCR
- Multiplex amplification
- High sensitivity applications
- Low-copy number PCR
- Complex DNA template PCR
- “Direct-blood” PCR
- “Anti-contamination” PCR
- “fast” PCR

General Protocol for amplification with HS-Fuzz DNA polymerase

The optimal reaction conditions for **HS-Fuzz** DNA polymerase differs from PCR protocols for standard (Taq-like) DNA polymerases. PCR conditions for **HS-Fuzz DNA Polymerase** is similar in PCR conditions to **Phusion (ThermoScientific)**, **Q5 polymerase (NEB)** or **iProof (Bio-Rad)** DNA polymerase, e.g. **HS-Fuzz** polymerase works better at elevated denaturation and annealing temperatures.

PCR reactions should be set up even at room temperature, because of blocking of enzymatic activity by **specific monoclonal antibodies** at PCR setup conditions (usually RT).

HS-Fuzz activates without pre-heating step, during "Initial Denaturation" PCR step. There is no need to enter additional heating step during PCR run.

Add and mix the following component on ice:

Component	50 μ L reactions	25 μ L reactions	Final concentration
PCR grade Water	Up to 50 μ L	Up to 25 μ L	
2.5x Uni Buffer-1500*	20 μ L	10 μ L	1X
2,5 mM MIX dNTPs	4 μ L	2 μ L	0.2 mM each
Primers			0.3-0.5 μ M each
Template DNA	optionally	optionally	
HS-Fuzz (2 U/ μ L)*	0,5-0,25 μ L	0.25-0,125 μ L	0.02U/ μ L

We recommend to use 50 μ l reaction for the PCR with HS-Fuzz polymerase

Reaction Components

1. HS-Fuzz Polymerase

An optimal amount of enzyme in 50 μ l reaction is 1U. In some cases it could be reduced up to 0.5U per reaction, depending on the length and complicity of amplified DNA sequence. For non-complex amplicons with the length less than 500bp and GC-content <60% amount of HS-Fuzz could be decreased up to 0.5-1.25U per 50 μ lreaction.

2. Buffer

2.5 X Uni HS-Fuzz buffers provides very high reproducibility across the wide range of amplification conditions, including "fast-PCR" (reduced time of PCR reaction). 2,5X Uni Buffer contains 1.5mM Mg²⁺, as the final concentration. In some cases we recommend to optimize Mg concentration in the range 1.5-4.5mM

3. DNTP's

For most of applications 200 μ M of each of dNTP's as final concentration is an optimal. It's not necessary to optimize dNTP's concentration. dUTP or other dUTP derivatives should be used replacing TTP in PCR reaction for "anti-contamination" PCR.

4. Primers

Usually 10-20pmol of each specific primer in reaction is enough to get good PCR result. If you are using 2-step PCR with the whole blood as a template, it's better to use \geq 20pmol of each primer.

5. PCR Additives

HS-Fuzz polymerase is compatible with the most of commonly used PCR additives for enhancing of high GC-content DNA templates (glycerol, betaine, DMSO and other).If one will use any additives; take into account the changes of Tm of primers and DNA to correct annealing temperature.

Cycling Protocol:

Cycle step	2-step amplification		3-step amplification		Cycles
	T $^{\circ}$ C	Time	T $^{\circ}$ C	Time	
Initial Denaturation	98$^{\circ}$C	1-5 min	98$^{\circ}$C	1-5min	1
Denaturation	98$^{\circ}$C	2-10 S	98$^{\circ}$C	2-10 S	25-35
Annealing	-	-	58-72*	10-30 S	
Extension	66-72$^{\circ}$C	15-30 S/Kb	72$^{\circ}$C	15-30 S/Kb**	
Final extension	72$^{\circ}$C	1-3 min	72$^{\circ}$C	1-3 min	1
	4$^{\circ}$C	hold	4$^{\circ}$C	hold	

*Optimal Tm for the primer pair recommended as Tm of the lower primer, for the standard oligos <20nt.

For **HS-Fuzz** polymerase Tm of the primers should be corrected, as +3-10 $^{\circ}$ C, comparing with "regular" Taq-based enzymes PCR conditions.

NOTE: If whole blood is using in PCR reaction as a template, it's highly recommended to hold primers concentration above 15 pmol (preferable 20-25 pmol) and to titrate Mg concentration from 2mM (1X) up to 3, 5-4mM, depending on the volume of blood in PCR.

1. Denaturation:

- 1) Initial denaturation for 5 min at 98°C is necessary only for blood cells lysis;
- 2) for most applications, including "direct-blood" 5 sec at 98°C is enough for denaturation of the sample in during PCR run

NOTE: Don't use lower T denaturation, then 98°C; it can cause problems in PCR (nonspecific amplification, poor yield of PCR product, etc.)

2. Annealing/Extension:

For **HS-Fuzz** polymerase "Annealing" and "Extension" steps should be combined if:

- T_m of both primers are not differs dramatically (<3°C);
- T_m of the primers are >65°C (optimal T_m for the primers lays between 65-70°C)

If primers T_m is about 60-61°C for both primers ones can apply simple formula to determine starting Ta/e point - (T_m of the lower primer +72°C)/2. For most applications it works fine.

To determine a better Ta/e run gradient amplification.

To avoid nonspecific band formation/smearing during amplification not exceed extension time of 30 seconds and use the highest ramp rate of amplificator (the ramp rate >4-5°C preferable)

**For non-complex DNA templates (plasmid DNA, phage DNA, BAC clone) extension time could be reduced up to 15 sec/Kb.

For complex DNA templates (human DNA) in amplification of the fragments longer than 2Kb strongly recommended to apply Extension time as 30 sec/Kb

For the DNA templates with GC content less 65% and up to 500bp length it's quite enough to set Extension time as 15 sec per cycle

3. Post-PCR amplicon detection:

- If you are using whole blood as a template, after amplification spine-down blood cells debris to avoid its application to the gel for DNA detection.

Note: Don't use excess of the blood (>10%), because the post-PCR debris volume of blood cells is very high, and not allow to run more than 1-2 gels.