

Data sheet

This data sheet provided for information only

SmarTaq DNA POLYMERASE (Low Glycerol)**RESEARCH USE ONLY****Recombinant SmarTaq DNA Polymerase**

(Deoxynucleosidetriphosphate: DNA Deoxynucleosidyltransferase E.C. 2.7.7.7.)

SOURCE:**Thermus Auqaticus YT1 strain, Mouse monoclonal antibodies.**

Cat. No.	Pack	Conc.
DSTL-500	500 Units	5 U/μl
DSTL-1000	1000 Units	5 U/μl

DESCRIPTION

SmarTaq DNA Polymerase is complex mixture of a thermostable 94 kDa **DiaTaq** DNA Polymerase purified from E.coli PVG-AI recombinant strain expressing **Thermus aquaticus** polymerase gene and **highly specific monoclonal antibodies**.

SmarTaq is inactive under conditions of amplification reaction preparation. It provides **improved specificity and high yield amplification, even with low copy number and complex DNA templates**, when compared to standard DNA polymerases. **SmarTaq** can **eliminate amplification artifacts**, such as primer-dimer formation and mispriming.

An advantage of **SmarTaq** is the absence of additional heating step for polymerase activation. Heat activation of enzyme occurs during the first denaturation step. An active complex of **SmarTaq** dissociates automatically over **+70°C**, allowing activating DNA polymerase, preventing mispriming during amplification.

Stability:

Shelf life 24 months if store at -20°C

UNIT DEFINITION

One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTP's into acid-insoluble form in 30 minutes at 74°C under assay conditions.

STORAGE AND DILUTION BUFFER:

20mM Tris-HCL (pH 8.0);100mM KCL;0.1mM EDTA; 1mM DTT; 5% glycerol, 0.5% Nonidet P-40;0.5% Tween-20, Stabilizers.

AMPLIFICATION BUFFERS:

10X NH4-buffer: 166mM (NH₄)₂SO₄; 670mM Tris-HCL (pH 8.8 at 25°C); 0.1% Tween-20.

5X "Ready-to-load" NH4-buffer: 166mM (NH₄)₂SO₄; 670mM Tris-HCL (pH 8.8 at 25°C);0.1% Tween-20; two inherent dyes; stabilizers

5X "Enhanced" NH4-buffer:

166mM (NH₄)₂SO₄; 670mM Tris-HCL (pH 8.8 at 25°C); 0.1% Tween-20; stabilizers/enhancer

ASSOCIATED ACTIVITIES:

Endonuclease and exonuclease activities were not detectable after 2 and 1 hours incubation, respectively, of 1 μg lambda DNA and 0.22 μg of EcoR I digested lambda DNA, respectively, at 72°C in the presence of 15-20 units of **SmarTaq** DNA polymerase

APPLICATIONS:

- Real-Time PCR
- Multiplex amplification - **Highly recommended**
- High sensitivity applications
- Low-copy number PCR
- Complex DNA template PCR
- Preparation of lyophilized enzyme or MasterMixs

STORAGE CONDITIONS :

Store **SmarTaq** DNA Polymerase at -20°C for long term storage

Can be stored for 3 months at +4°C

SHIPPING CONDITIONS:

Should be shipped at ambient temperature

For long distance shipments preferably in

Blue Ice